Istruzioni d'uso

ver. 3 - 20/04/2022

Biologia Molecolare

DPYD Real Time (FRET)

REF

AA1579/48A

48 TEST

CND

W0106010499

CE

IVD

NUCLEAR LASER MEDICINE S.r.l.

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
IVD

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The **AA1579/48A** device provides the material for the analysis in DNA extracted from human blood samples of five DPYD (dihydropyrimidine dehydrogenase) gene variants that reduce enzymatic activity of the DPD protein involved in the catabolism of fluoropyrimidines:

c.1905+1G>A, c.1679T>G, c.2846A>T, c.1129–5923C>G, c.2194G>A.

The device is based on amplification of the target sequences and analysis of melting curves.

A mix for each variant, including amplification enzyme and Uracil-DNA Glycosylase (UNG) to prevent amplified contamination, is provided.

This device has been validated with:

- manual column based kit (NLM code AA1001), extraction kit for automatic instruments OMNIA LH75 (NLM code AA1319/48) and with the most common automatic nucleic acid extraction systems.
- CFX Real Time PCR (BioRad) instrument.

INTRODUCTION

Fluoropyrimidines (5-fluorouracil, capecitabine and tegafur) are a drugs class used in the treatment of many types of solid tumors. Metabolism of fluoropyrimidines requires DPD protein (dihydropyrimidine dehydrogenase) encoded by DPYD gene. DPD catabolyses 85% of the administered dose of fluoropyrimidines and reduced or absent activity of this enzyme can result in severe and sometimes fatal toxicity.

Functional DPYD gene variants have been found to be associated with reduced or abrogate DPD activity.

Recommendations for pharmacogenetic analysis of DPD and fluoropyrimidines are presented in 2019 from AIOM-SIF in which is recommended a test to analyze 5 mutations with clinical relevance to all patients candidates for treatment with fluoropyrimidines:

| HGVS | SNP | Allele | Effetto |
|------------------|------------|---------|------------|
| c.1905+1G>A | rs3918290 | DPYD*2A | IVS14+1G>A |
| c.1679T>G | rs55886062 | DPYD*13 | p.I560S |
| c.2846A>T | rs67376798 | | p.D949V |
| c.1129–5923C>G ^ | rs75017182 | | IVS10 C>G |
| c.2194G>A | rs1801160 | DPYD*6 | p.V732I |

^ Variant c.1129–5923C>G is in perfect linkage disequilibrium with c.1236G> A (rs56038477) in European population.

Fluoropyrimidines dose omission or reduction were recommended in patients carriers of these variants:

| DPYD Genotype | | Recommended fluoropyrimidines dose |
|---------------|------------|------------------------------------|
| Wild Type | c.1236GG | 100% |
| | c.1679TT | |
| | c.1905+1GG | |
| | c.2846AA | |
| | c.2194GG | |
| Heterozygous | c.1236GA | 75% |
| | c.1679TG | 50 |
| | c.1905+1GA | |
| | c.2846AT | |
| | c.2194GA | 85% |
| Homozygous | c.1236AA | 50% |
| | c.1679GG | No fluoropyrimidines |
| | c.1905+1AA | |
| | c.2846TT | |
| | c.2194AA | 70% |


REFERENCES

- Amstutz *et al*, Clin Pharmacol Ther 2018; 103(2): 210–216
- AIOM-SIF, Raccomandazioni 2019 per analisi farmaco genetiche
- Boige *et al*, JAMA Oncol 2016; 2: 655-662
- Del Re *et al*, Pharmacogenomics J 2019; 19: 556-563
- Henricks *et al*, Lancet Oncol 2018;19:1459-1467
- Iachetta, Bonelli *et al*, British Journal of Cancer 2019; 120: 834-839
- Meulendijks *et al*, Cancer Chemotherapy and Pharmacology 2016; 78: 875–880
- Offer *et al*, Cancer Res 2013; 73(6): 1958-1968

PRINCIPLES OF THE PROCEDURE

The assay is based on DNA amplification and Fluorescent Resonance Energy Transfer (FRET) that occurs between two adjacent probes complementary to the region that includes the mutation: a 5' probe includes the mutation of interest and it is labeled with a donor fluorophore; instead a 3' probe carries the acceptor fluorophore. If the two probes appear in their target sequence, the two fluorophores come to be in proximity and, in the presence of a light source, energy is therefore transferred from the donor to the acceptor; in the absence of the target, energy transfer does not take place.

The probes allow to identify variants of the target sequence by means of the dissociation curve (melting curve); the amplicon is subjected to a rising temperature ramp with consequent detachment of the probes from the DNA. If the 5' probe is complementary to the mutated DNA strand, the hybridization between the two sequence will be more stable than the duplex between the probe itself and the strand that is not mutated (single point mismatch), with consequent peak of dissociation to an higher temperature. Therefore, through the melting curve analysis it is possible to distinguish the three genotypes based on the peak positions: two single peak with different temperatures for the homozygous samples and double peak for the heterozygote.

PRODUCT COMPOSITION

Store at -25/-15°C

| Reagents | NLM code | N° vial | Volume/vial |
|--|----------|---------|-------------|
| Mastermix DPYD*2A , red cap Solution containing oligonucleotides, fluorescent probes and enzymes | KA1278-n | 1 | 600 µl |
| Mastermix DPYD*13 , green cap Solution containing oligonucleotides, fluorescent probes and enzymes | KA1279-n | 1 | 600 µl |
| Mastermix DPYD D949V , blue cap Solution containing oligonucleotides, fluorescent probes and enzymes | KA1280-n | 1 | 600 µl |
| Mastermix DPYD IVS10 , yellow cap Solution containing oligonucleotides, fluorescent probes and enzymes | KA1281-n | 1 | 600 µl |
| Mastermix DPYD V732I , purple cap Solution containing oligonucleotides, fluorescent probes and enzymes | KA1282-n | 1 | 600 µl |
| DNase, RNase, Protease-free water , transparent cap | AZ040-n | 1 | 1600 µl |

The amount of reagents indicated in the table is referred to the standard size of the kit.
Reduced packaging is available on request for evaluation and/or demonstration of the product.

STABILITY AND STORAGE

- All the reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (-25/-15°C)
- Thaw Mastermix on ice or at +2/+8°C
- Probes in the Mastermix are photosensitive: avoid prolonged exposure to light
- At the end of each assay, store reagents at their correct storage temperature.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established Good Laboratory Practices and universal precautions; wear personal protective apparel
- Discard the used disposable material, worn gloves and all reagents as special waste
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled
- If eyes, skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice
- Don't use components beyond the expiration date
- Don't mix reagents from different lots
- Keep reagents separated to avoid possible contamination (samples and amplification products)
- Avoid repeated freezing and thawing of the mixes
- It is recommended to make analysis in three different areas:
 - Area 1: samples handling and extraction; automatic PCR setup
 - Area 2: manual PCR setup
 - Area 3: post PCR (Real Time PCR)
- Don't use the device if the box is damaged; contact the supplier
- It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct work.

MATERIALS REQUIRED BUT NOT SUPPLIED

- Automatic DNA Extractor and disposables
- Workstation for PCR set-up and disposables
- Vertical downflow airbox
- DNA extraction kit
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR strips or plates
- Real Time PCR instruments
- Heterozygous positive control, single use (cod. NLM FA157/10).

For compatibility with instruments other than that indicated, contact *Nuclear Laser Medicine srl*.

PROCEDURE

Use fresh or correctly stored +2/+8°C (for up to 2 days) or -25/-15°C for longer period blood samples. Use only EDTA or citrate as anticoagulants (do not use heparin).

DNA EXTRACTION

- Manual isolation kit (NLM code AA1001)
- extraction kit (NLM code AA1319/48) for OMNIA LH75 automatic system and the most common automatic nucleic acid extraction systems.

For the use of OMNIA LH75 automatic system, follow the instructions provided by *Nuclear Laser Medicine srl*; for the extraction with other automatic systems, follow the supplier indications.

AMPLIFICATION

Warning

- Avoid using blood samples after repeated freeze/thaw cycles or not correctly stored in order to prevent bad assay results
- Before use, gently shake the Mastermix and check the absence of residues in the cap
- Extracted DNA concentration must be about 15-60 ng/μl.
- The Positive Control is for single use only. Discard the remaining material in the tube.

manual PCR setup

For each of the five Mastermixes, proceed as follows:

| Reagents | Volume/sample |
|-----------------------------------|---------------|
| Mastermix DPYD | 10 μl |
| DNase, RNase, Protease-free water | 6 μl |
| DNA | 4 μl |

Mix well the three reagents together

If possible, spin strips or plate before amplification to eliminate the presence of bubbles.

PCR setup with automatic system

For PCR set up with the automatic OMNIA LH75 system, use the protocol provided by *Nuclear Laser Medicine srl*.

THERMAL PROFILE SETTING AND RESULTS ANALISYS

 For **CFX96 (Biorad)** instrument use, download the Instructions for Use (NLM code **DX007**) from the www.nlm.it. For compatibility with other Real Time PCR instruments, contact *Nuclear Laser Medicine srl*.

THERMAL PROFILE

Setting up *Sample Volume*: 20 µl
Enter the following thermal profile:

- 1 45,0° C for 5:00
- 2 95,0° C for 2:00
- 3 95,0° C for 0:15
- 4 54,0° C for 0:20 + Plate Read
- 5 72,0° C for 0:10
- 6 **GOTO 3, 39 more times**
- 7 35,0° C for 2:00
- 8 75,0° C for 1:00, Slow Ramp Rate to 0,1°C per second
- 9 25,0° C for 10:00
- 10 Melt Curve 35,0 to 75,0° C, increment 0,5° C for 0:05 + Plate Read


Figure 1. Reproduction of the thermal profile, as displayed in the instrument program.

Plate Setup

- Select *Sample Type* → *Unknown* for the samples to be run
- Select *Fluorophores* → FAM
- Select *Target Name* and enter the description shown in the table below:

| Target | Fluorophore |
|------------|-------------|
| DPYD*2A | FAM |
| DPYD*13 | |
| DPYD D949V | |
| DPYD IVS10 | |
| DPYD V732I | |

- For Positive controls, select: - *Sample Type* → *Positive Control*
- *Sample Name* → *ctrl DPYD etero*
- For Negative controls (No Template Control), select *Sample Type* → *NTC* (*Sample name* of your choosing)

INTERPRETATION OF RESULTS

At the end of the run, the *Data Analysis* window opens automatically.

- Select the *Melt Curve* window: the graph on the left shows the raw data of the melting curve while the graph on the right shows the peaks of the melting curves (analysis -d (RFU)/dT)
- Set up *Peak Type* → *Negative*

The instrument automatically selects the threshold (blue line): move it at will to eliminate background noise and/or non-specific signals

For each sample, analyze the melting peaks selecting the individual wells on the plate.

Warning: Do not consider peaks out of the temperature range.

| Target | T _m | |
|------------|----------------|------------|
| | Wild tipe | Mutant |
| DPYD*2A | 50.5°C ± 2 | 57.0°C ± 2 |
| DPYD*13 | 59.0°C ± 2 | 52.0°C ± 2 |
| DPYD D949V | 52.0°C ± 2 | 57.0°C ± 2 |
| DPYD IVS10 | 48.5°C ± 2 | 57.5°C ± 2 |
| DPYD V732I | 48.0°C ± 2 | 58.5°C ± 2 |


If both peaks are present (wild type and mutant), sample should be considered heterozygous.

Important: if heterozygous samples are found, evaluate the raw data in order to exclude the inclusion of non-specific peaks.


⇒ For automatic interpretation of results and report generation, use the **Real-PHARMA** software following the NLM code **DO041** user manual or RealGene software, NLM **DO022** code.

The Positive Control inclusion in the run is optional. Non-compliant CP results do not invalidate the run.


Warning: the amplification curves do not contain useful information for genotyping which is based only on peaks T_m.


Wild type and heterozygous samples analyzed with **DPYD*2A** mix


Wild type and heterozygous samples analyzed with **DPYD*13** mix


Wild type and heterozygous samples analyzed with **DPYD D949V** mix

Wild type and heterozygous samples analyzed with **DPYD IVS10** mix

Wild type and heterozygous samples analyzed with **DPYD V732I** mix

Wild type sample for all five variants and threshold selection

Figure 2. Wild type and heterozigous samples analyzed with DPYD Mastermixes

WARNING

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of genetic mutation in the primer and/or probe specific target regions.

Figure 3 share the T_m of rs17376848 and rs3918289 SNPs that are close to DPYD*2A mutation and which have a relevant frequency in the population:


Figure 3. Shift of wild type peak (red) in presence of SNPs rs17376848 (green) e rs3918289 (blue). Result obtained using synthetic targets.


Figure 4. (A) Heterozygous sample for SNPs rs17376848 (red) compared with wild type (blue); (B) synthetic samples heterozygous for SNPs rs3918289A and (C) rs3918289C (red) compared with wild type (blue).

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

The analytical sensitivity of "DPYD Real time (FRET)" kit, expressed as the minimum quantity of the target that can be detected, is 5 ng/µl of DNA (equivalent to 20 ng/reaction).

Diagnostic Sensitivity

Diagnostic sensitivity of "DPYD Real Time (FRET)" kit was assessed by analyzing 30 DNA samples with different genotypes.

All samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Diagnostic Specificity

The specificity "DPYD Real Time (FRET)" kit was determined by analyzing 30 DNA samples with different genotypes. No false positive results were obtained.

In agreement with these data the specificity of the kit is 100%.

Reproducibility

Intra-assay reproducibility

Intra-assay reproducibility was assessed by analyzing 3 samples of different genotype in 4 replicates in the same run and on different instruments.

Inter-assay reproducibility

Inter-assay reproducibility was assessed by analyzing 3 samples of different genotype in 4 replicates in different run and on different instruments.

On the basis of the results obtained the reproducibility is 100%.