Istruzioni d'uso

ver. 2 - 19/04/2022

Biologia Molecolare

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UGT1A1 Real Time (FRET)

REF

AA1628/24A

CE

24 TEST

CND

W010699

IVD

NUCLEAR LASER MEDICINE S.r.l.

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
60126

NUCLEAR LASER MEDICINE S.r.l.

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY
Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: www.nlm.it - E-MAIL: segreteria@nlm.it

ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).

INTENDED USE

The **AA1628/24A** device provides the material for the analysis in human genomic DNA samples of four allelic variants present in the UGT1A1 (UDP-glucuronosyltransferase 1A1 isoform) gene promoter: UGT1A1*1, UGT1A1*28, UGT1A1*36 e UGT1A1*37.

The device is useful before to start a therapy with irinotecan or to confirm the diagnosis of Gilbert's syndrome.

The device is composed by one mix (including amplification enzyme) and it is based on amplification of the target sequences and analysis of melting curves.

This device has been validated with:

- manual column based kit (NLM code AA1001), extraction kit for automatic instruments OMNIA LH75 (NLM code AA1319/48) and with the most common automatic nucleic acid extraction systems
- CFX Real Time PCR (BioRad) instrument.

INTRODUCTION

Irinotecan (brand name Camptosar) is a drug widely used in the treatment of advanced or metastatic colorectal cancer. Irinotecan therapy is associated with a high incidence of toxicity, including severe neutropenia and diarrhea. Irinotecan is inactivated by a UDP-glucuronosyltransferase (UGT) enzyme encoded by the UGT1A1 gene that transforms lipophilic metabolites (such as bilirubin, steroids, toxins, and drugs) into water-soluble metabolites that can be excreted from the body.

UGT1A1 gene transcriptional activity is affected by the length of the TATA box in the promoter, because a larger number of TA repeats results in reduced transcriptional activity:

- UGT1A1*1** [(TA)₆TAA] is the wild type allele associated with normal enzyme activity
- UGT1A1*28** [(TA)₇TAA] variant contains an extra thymine-adenine (TA) repeat associated with a decrease of transcriptional activity of the gene; the presence of two copies of this variant results in about a 70% decrease that is associated with Gilbert's syndrome (a mild form of conjugated hyperbilirubinemia) and predisposition to severe toxicity of irinotecan
- UGT1A1*37** [(TA)₈TAA] characterized by 8 TA repetitions, it is associated with a reduced enzymatic activity comparable to that observed in the presence of the allele * 28
- UGT1A1*36** [(TA)₅TAA] allele, has only 5 TA repeats and is associated with increased promoter activity with reduced risk to develop irinotecan toxicity.

*36 e *37 variants occur almost exclusively in African population while *28 allele is commonly found in African-American and Caucasian population.

AIOM-SIF document, presented in 2019, recommend pharmacogenetic analysis of UGT1A1 in order to prevent potentially serious adverse reactions to irinotecan therapy. If the homozygous presence of the variant * 28 is detected, a reduction in the dose of drug administered is in fact recommended:

UGT1A1 Genotype		Recommended irinotecan dose
Wild-type	*1/*1	100%
Heterozygous	*1/*28	100%
Mutated homozygous	*28/*28	70%

For the customization of irinotecan therapy, the allele *36 can be interpreted as an *1 allele and allele *37 as an allele *28.


REFERENCES

- AIOM-SIF, Raccomandazioni 2019 per analisi farmaco genetiche
- Barbarino *et al*, Pharmacogenet Genomics 2014; 24(3): 177–183
- Cecchin *et al*, Drug Resist Updat 2018; 39:18-40
- Deal, Medical Genetics Summaries, updated 2018; 1-12
- Shulman *et al*, Cancer 2011; 117(14): 3156–3162
- Takano and Sugiyama, Pharmacogenomics and Personalized Medicine 2017; 10: 61–68
- <https://www.pharmgkb.org/guideline/PA166127626>.

PRINCIPLES OF THE PROCEDURE

The assay is based on DNA amplification and Fluorescent Resonance Energy Transfer (FRET) that occurs between two adjacent probes complementary to the region that includes the mutation: a 5' probe includes the mutation of interest and it is labeled with a donor fluorophore; instead a 3' probe carries the acceptor fluorophore. If the two probes appear in their target sequence, the two fluorophores come to be in proximity and, in the presence of a light source, energy is therefore transferred from the donor to the acceptor; in the absence of the target, energy transfer does not take place.

The probes allow to identify variants of the target sequence by means of the dissociation curve (melting curve); the amplicon is subjected to a rising temperature ramp with consequent detachment of the probes from the DNA. The 5' probe is complementary to the [(TA)₈TAA] variant; the hybridization between the two sequences will be more stable than the duplex between the probe itself and the strand with [(TA)₇TAA], [(TA)₆TAA] or [(TA)₅TAA] variants, with consequent peak of dissociation to an higher temperature. Therefore, through the melting curve analysis it is possible to distinguish the four different genotypes based on the peak positions: a single peak with different temperatures for the homozygous samples and double peak for the heterozygote.

PRODUCT COMPOSITION

Store at -25/-15°C

Reagents	NLM Code	N° vial	Volume/vial
MasterMix UGT1A1 , red cap Solution containing oligonucleotides, fluorescent probes and enzymes	KA1313-n	1	470 µl
DNase, RNase, Protease-free water , transparent cap	AZ040-n	1	150 µl

The amount of reagents indicated in the table is referred to the standard size of the kit.

Reduced packaging is available on request for evaluation and/or demonstration of the product.

STABILITY AND STORAGE

- All the reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (-25/-15°C)
- Thaw MasterMix on ice or at +2/+8°C
- Probes in the MasterMix are photosensitive: avoid prolonged exposure to light
- At the end of each assay, store reagents at their correct storage temperature.

PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit. Handle this product according to established Good Laboratory Practices and universal precautions; wear personal protective apparel
- Discard the used disposable material, worn gloves and all reagents as special waste
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled
- If eyes, skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice
- Don't use components beyond the expiration date
- Don't mix reagents from different lots
- Keep reagents separated to avoid possible contamination (samples and amplification products)
- Avoid repeated freezing and thawing of the mixes
- It is recommended to make analysis in three different areas:
 - Area 1: samples handling and extraction; automatic PCR setup
 - Area 2: manual PCR setup
 - Area 3: post PCR (Real Time PCR)
- Don't use the device if the box is damaged; contact the supplier
- It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct work.

MATERIALS REQUIRED BUT NOT SUPPLIED

- Automatic DNA Extractor and disposables
- Workstation for PCR set-up and disposables
- Vertical downflow airbox
- DNA extraction kit
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR strips or plates, clear
- Real Time PCR instruments
- Heterozygous positive control, single use (cod. NLM FA165/10).

PROCEDURE

Use fresh or correctly stored +2/+8°C (for up to 2 days) or -25/-15°C for longer period blood samples. Use only EDTA or citrate as anticoagulants (do not use heparin).

DNA EXTRACTION

- Manual isolation kit (NLM code AA1001)
- Extraction kit (NLM code AA1319/48) for OMNIA LH75 automatic system and the most common automatic nucleic acid extraction systems.

For the use of OMNIA LH75 automatic system, follow the instructions provided by *Nuclear Laser Medicine srl*; for the extraction with other automatic systems, follow the supplier indications.

AMPLIFICATION

Warning

- Avoid using blood samples after repeated freeze/thaw cycles or not correctly stored in order to prevent bad assay results
- Before use, gently shake the MasterMix and check the absence of residues in the cap
- Extracted DNA concentration must be about 15-60 ng/μl.
- The Positive Control is for single use only. Discard the remaining material in the tube.

Manual PCR setup

For the manual PCR set up, proceed as follows:

Reagents	Volume/sample
MasterMix UGT1A1	13 μl
DNase, RNase, Protease-free water	3 μl
DNA	4 μl

Mix well the three reagents together.

If possible, spin strips or plate before amplification to eliminate the presence of bubbles.

PCR setup with automatic system

For PCR set up with the automatic OMNIA LH75 system, use the protocol provided by *Nuclear Laser Medicine srl*.

THERMAL PROFILE SETTING AND RESULTS ANALYSIS

CFX Real Time PCR (Biorad)

For **CFX96 (Biorad)** instrument use, download the Instructions for Use (NLM code **DX007**) from the www.nlm.it.

THERMAL PROFILE

Setting up *Sample Volume*: 20 µl

Enter the following thermal profile:

- 1 95,0° C for 2:00
- 2 95,0° C for 0:20
- 3 60,0° C for 0:20 + Plate Read
- 4 72,0° C for 0:10
- 5 **GOTO 2, 49 more times**
- 6 95,0° C for 1:00
- 7 30,0° C for 3:00
- 8 Melt Curve 30,0 to 60,0° C, increment 0,3° C for 0:05 + Plate Read


Figure 1. Reproduction of the thermal profile, as displayed in the instrument program.

Plate Setup

- Scan mode → SYBR/FAM only
- Select Fluorophores → FAM
- Select Target Name and enter the description shown in the table below:

Target	Fluoroforo
UGT1A1	FAM

- Select Sample Type → Unknown for the samples to be run
- Select Sample Type → Positive Control for Positive controls
Sample Name → ctrl1-28 or ctrl36-37
- Select Sample Type → NTC for Negative controls (No Template Control)

INTERPRETATION OF RESULTS

At the end of the run, the *Data Analysis* window opens automatically.

- Select the *Melt Curve* window: the graph on the left shows the raw data of the melting curve while the graph on the right shows the peaks of the melting curves (analysis -d (RFU)/dT)
- Set up *Peak Type* → Negative

The instrument automatically selects the threshold (blue line): move it at will to eliminate background noise and/or non-specific signals

For each sample, analyze the melting peaks selecting the individual wells on the plate.

Target		Tm
UGT1A1*36	[(TA) ₅ TAA]	39°C ±1.5
UGT1A1*1	[(TA) ₆ TAA]	43,5°C ±1.5
UGT1A1*28	[(TA) ₇ TAA]	47°C ±1.5
UGT1A1*37	[(TA) ₈ TAA]	51,5°C ±1.5

If two peaks are present, sample should be considered heterozygous.

Important: if heterozygous samples are found, evaluate the raw data in order to exclude the inclusion of non-specific peaks.


For automatic interpretation of results and report generation, use the **Real-PHARMA** program following the NLM code **DO041** user manual.


Positive Control inclusion in the run is optional. Non-compliant CP results do not invalidate the run.

Warning

- The amplification curves do not contain useful information for genotyping which is based only on peaks Tm
- Do not consider peaks out of the temperature range
- As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of genetic mutation in the primer and/or probe specific target regions.

Peaks of four UGT variant: *1, *28, *36 and *37

Homozygous sample for UGT1A1*1 e UGT1A1*28 (blue); heterozygous sample UGT1A1*1/*28 (red)

Figure 2. Homozygous and heterozygous samples analyzed with UGT1A1 MasterMix

TROUBLESHOOTING

Heterozygous sample UGT1A1*1/*28 and UGT1A1*1/*36

It's possible that heterozygous samples *1/*28 or *1/*36 share an unbalanced double peak (red sample, Figure 3). This result could be associated with blood or DNA sample of low quality. It is recommended to use sample thawed only few times and to evaluate raw data before to assign the genotype.


Figure 3. UGT1A1*1 and UGT1A1*28 samples (blue); *1/*28 sample (red)

PERFORMANCE CHARACTERISTICS

Analytical Sensitivity

The analytical sensitivity of "UGT1A1 Real Time (FRET)" kit, expressed as the minimum quantity of the target that can be detected, is 1 ng/μl of DNA (equivalent to 4 ng/reaction).

Diagnostic Sensitivity

Diagnostic sensitivity of "UGT1A1 Real Time (FRET)" kit was assessed by analyzing 36 DNA samples with different genotypes.

All samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

Diagnostic Specificity

The specificity "UGT1A1 Real Time (FRET)" kit was determined by analyzing 36 DNA samples with different genotypes. No false positive results were obtained.

In agreement with these data the specificity of the kit is 100%.

Reproducibility

Intra-assay reproducibility

Intra-assay reproducibility was assessed by analyzing samples with different genotype in 3 replicates in the same run.

Inter-assay reproducibility

Inter-assay reproducibility was assessed by analyzing samples with different genotype in 3 different runs.

On the basis of the results obtained the reproducibility is 100%.

LOT

Lot Number

REF

Catalogue
number

Temperature
Limitation

Use by

Keep away
from sunlight

Consult Instructions
for use

European
Conformity

Manufacturer

IVD


in vitro diagnostic
device

Sufficient for

CONTROL +

Positive
Control

Modified text
compared with
the previous
version