



Papilloplex HR-HPV mRNA Kit

Catalogue # (REF) : MPAHPV004

FOR *IN VITRO* DIAGNOSTIC USE

Store at -20°C

Protect from light

Instructions for Use – English

Version 4

Table of Contents

| | | |
|-------|---|----|
| 1 | Trademarks | 3 |
| 2 | Copyright..... | 3 |
| 3 | Kit contents | 3 |
| 4 | Shipment and storage | 3 |
| 5 | Introduction | 4 |
| 5.1 | Intended use | 4 |
| 5.2 | Target environment | 4 |
| 5.3 | Principle | 4 |
| 5.4 | Positive Controls provided in the kit | 5 |
| 6 | Contraindications, warnings and precautions | 5 |
| 6.1 | Contraindications | 5 |
| 6.2 | Interfering substances | 5 |
| 6.3 | Warning and precautions | 5 |
| 7 | Operating procedure..... | 6 |
| 7.1 | Specimen collection | 6 |
| 7.2 | RNA/DNA extraction | 6 |
| 7.3 | PCR Reaction Mix setup | 6 |
| 7.4 | Real-Time PCR instrument settings | 7 |
| 7.4.1 | Real-Time PCR instrument settings for Bio-Rad CFX96 | 7 |
| 7.4.2 | Real-Time PCR instrument settings for SLAN96P | 12 |
| 8 | Troubleshooting | 18 |
| 9 | Performance..... | 19 |
| 9.1 | Stability | 19 |
| 9.2 | Limit of Detection (LOD) | 19 |
| 9.3 | Analytical Specificity (cross-reactivity with others HPV types) | 20 |
| 9.4 | Analytical Specificity (cross-reactivity with other microorganisms) | 20 |
| 9.5 | Precision | 21 |
| 9.5.1 | Precision Intra-Assay | 21 |
| 9.5.2 | Precision Inter-Assay | 22 |
| 9.6 | Clinical Performance | 24 |
| 9.6.1 | Clinical specimens..... | 24 |
| 9.6.2 | Positive rate of HPV E6/E7 mRNA assay by histological diagnosis..... | 24 |
| 9.6.3 | Concordance between Papilloplex® HR-HPV mRNA kit and Cytology diagnosis..... | 25 |
| 9.6.4 | Correlation of Papilloplex® HR-HPV mRNA kit or Cytology with Histological diagnoses (CIN2+) | 25 |
| 10 | Product Specification | 26 |
| 11 | Symbols | 26 |
| 12 | Customer contact information..... | 27 |

1 Trademarks

Papilloplex® is a registered mark of GeneFirst Ltd registered in the UK and Ireland.

ThinPrep® is a registered mark property of Hologic Inc.

DNA AWAY™ is a trademark property of Molecular Bio-Products Inc.

2 Copyright

This document is property of GeneFirst Ltd including without limitation, all text, formats, graphics and logos and are protected from unauthorized copying and dissemination by the Copyright, Designs and Patents Act 1988 (as amended), by various intellectual property laws and by international conventions.

3 Kit contents

Materials supplied with the kit:

| Tube colour | Tube cap colour | Reagent | Description |
|-------------|-----------------|------------------|----------------------------------|
| Transparent | Blue | RT Mix | Reverse Transcriptase enzyme |
| Transparent | Green | Enzyme Mix | Taq Polymerase, buffer and dNTPs |
| Amber | Amber | Working Mix | Primers and probes |
| Transparent | Orange | Positive control | Control DNA |

Additional equipment & reagents required (not provided in the kit):

- Reagents and equipment for specimen collection, filtration, and DNA extraction
- Water, distilled (molecular biology grade)
- DNase, RNase and human DNA-free pipette tips with aerosol barriers
- DNase, RNase and human DNA-free tubes for preparing Reaction Mix
- Pipettes (adjustable)
- Tube racks
- Vortex mixer
- Microcentrifuge
- Real-Time PCR System: Clinical Performance evaluation have been performed on Bio-Rad CFX96 Real-Time PCR detection Systems. Analytical performance evaluation has been performed on both Bio-Rad CFX96 Real-Time PCR detection Systems and Shanghai Hongshi MedicalSlan-96P Real-Time PCR System.
- Real-Time PCR System Sequence Detection Software: CFX Manager software-IVD v1.6 and SLAN-96P 8.2.2
- PCR tubes, plates and accessories compatible with the use of the Real-Time PCR System
- Disposable powder-free gloves and lab-coat.

4 Shipment and storage

- Papilloplex® HR-HPV mRNA kit is shipped using frozen gel packs.
- Upon receipt of the kit, components must be stored in freezer at -20°C or below.
- The contents must be protected from light to prevent photobleaching and stored in the manufacturer's packaging.
- After opening, the kit is stable up to the expiration date indicated on the packaging provided that the components have been stored correctly according to the recommendations.

5 Introduction

Human papillomavirus (HPV) is one of the most common sexually transmitted infections and high-risk (HR) types of HPV cause the majority of cervical cancer cases. Based on their frequency in cervical cancer, 14 types (types 16, 18, 31, 33, 35, 39, 45, 51, 52, 56, 58, 59, 66 and 68) has been defined as carcinogens and, therefore, named as high-risk HPV (HR-HPV) types. The viral *E6* and *E7* oncogenes are active in cervical carcinomas and their corresponding proteins are directly involved in triggering cell proliferation, inhibition of apoptosis, reprogramming of differentiation, and chromosomal instability leading to malignant transformation of host cells. Transcripts of these oncogenes have always been observed in cervical carcinomas and are used in this assay as a predictive biomarker for the progression of HPV infection to cervical cancer. This kit complements the existing Papilloplex® HR-HPV Kit and Papilloplex® HR-HPV DNA Kit.

GeneFirst has developed the Multiplex Probe Amplification (MPA) technology enabling real-time PCR detection of multiple targets in a single closed-tube reaction. The Papilloplex® HR-HPV mRNA kit can detect and differentiate E6/E7 mRNA of all 14 HR HPV types in addition to a cellular control target in a single reaction.

5.1 Intended use

The Papilloplex® HR-HPV mRNA kit is a single reaction qualitative *in vitro* multiplex real-time polymerase chain reaction-based test for the screening of viral induced cellular oncogenic transcripts E6/E7 mRNA of 14 individual high-risk (HR) Human Papillomavirus (HPV) types in human cervical samples. For professional use only.

5.2 Target environment

The Kit is for professional use only and may be used in pathology laboratories or in epidemiological studies in research settings.

5.3 Principle

The Papilloplex® HR-HPV mRNA kit is based on the MPA technology patented by GeneFirst Limited. The MPA technology allows differentiation of up to six different targets per fluorescence channel, using a combination of PCR primers and probes (dual labelled fluorescent probe and partially complementary oligo hybrid) for each specific HR HPV target. Each probe has a unique melting profile (different melting temperature) that allows specific detection of the target present in the sample (Figure 1).




Figure 1. Melting curves/profiles in FAM channel. Y axis = derivative fluorescence and X axis = temperature. Red line denotes Negative Control - (melting profile of No Template Control, NTC) and green line denotes test sample/specimen. The difference (shown by grey circle) in the melting curve allows identification of the HR HPV genotype present in the sample tested (HR HPV type 16 in the figure).

Real-time reverse transcription PCR is performed in a single step thereby converting into cDNA and amplifying the E6/E7 mRNA of any HR HPV targets present in the sample. If mRNA of one or more targets are present, the corresponding probes are consumed during the amplification, producing a reduction in the signal visible in the

final melt curve analysis. Comparing melting profiles of the probes reveals which specific probe is consumed, thus indicating which HR HPV target is present in a sample.

5.4 Positive Controls provided in the kit

The tube with an **orange** cap contains Positive control to be used in each PCR run. The control contains HR HPV genotypes 16, 45, 18 and HPRT1 synthetic. They are detected in the FAM, JOE (HEX), ROX and Cy5 channels respectively.

6 Contraindications, warning and precautions

6.1 Contraindications

There are no known contraindications identified.

6.2 Interfering substances

The performance of the kit may be adversely affected by known PCR inhibitors co-extracted from patients' samples (such as blood, acetic acid, iodine, excessive mucous or pharmaceutical preparations such as lubricant gels, spermicide creams etc). The use of samples containing such substances should be avoided.

6.3 Warning and precautions

- This kit is designed to be used for *in vitro* diagnostic use and should be used by trained personnel with good laboratory practice and good competency in real-time PCR.
- Upon arrival, please check the kit for signs of damage. If damaged, please contact GeneFirst customer service or your local distributor. Do not use damaged kit components as they may not yield the expected performance.
- Partially used kit can be re-used if stored at -20°C following initial use. Thawing and freezing more than 4 times is not recommended. Prepared reaction mixes are for single use only and are not intended to be re-used.
- Do not use the product beyond its expiry date.
- Do not mix reagents from different batches.
- The positive control provided in the kit must be used as control in every experiments.
- Fluorescently labelled probes included in the Amber tube (Working Mix) are sensitive to photobleaching. Exposure to light should be avoided as much as possible.
- Lab coats and powder-free gloves must always be worn.
- Never touch the inside of the tube cap.
- Appropriate pipette tips with an aerosol barrier and free of DNase, RNase and human DNA must be used.
- Use appropriate measures to decontaminate working surfaces such as wipe/spray with 0.5% Sodium Hypochlorite solution or DNA AWAY™.
- Thaw all components thoroughly at room temperature before using the kit and mix.
- Avoid excessive vortexing.
- Disposal of unused reagents and waste must be done in accordance with country or local regulations.
- Safety Data Sheet (SDS) is available on request from either GeneFirst or your distributor.

7 Operating procedure

7.1 Specimen collection

The kit has been validated on LBC samples conserved in PreservCyt transport medium (ThinPrep®). For details on specimen collection, please refer to relevant product details from your supplier.

7.2 RNA/DNA extraction

Two DNA/RNA extraction kits have been tested for use in conjunction with the Papilloplex® HR-HPV mRNA kit Kit using their standard manufacturer recommended protocol:

- Quick DNA/RNA Viral Kit from Zymo Research
- Viral-PrepAdem-Kit from Ademtech

For detailed protocols on using individual kits, please refer to product details from your relevant suppliers.

7.3 PCR Reaction Mix setup

PCR Reaction Mix is prepared according to the table below. All steps can be performed at room temperature with minimal exposure to light. Before use, the Reaction Mix should be fully thawed and mixed thoroughly by vortexing and briefly spun down.

Transfer 15 µl of the Reaction Mix into each of the wells of a PCR plate, followed by adding 5 µl of positive control (PC) or test sample into each well. At least one PC should be included per run.

Note: adding consistent and precise amounts of reagents and DNA or control is critically important for accurate genotyping results.

| Tube Cap colour | Name | Volume per single reaction (µl) | Volume required for 96 reactions plus excess** (µl) |
|------------------------------------|--|---------------------------------|---|
| <i>Blue</i> | RT Mix | 1.00 | 100.00 |
| <i>Green</i> | Enzyme Mix | 4.00 | 400.00 |
| <i>Amber</i> | Working Mix | 2.00 | 200.00 |
| (Not supplied with the kit) | H ₂ O (molecular biology grade) | 8.00 | 800.00 |
| | Reaction Mix | 15 | 1500.00 |
| <i>And</i> | | | |
| <i>Orange*</i> | Positive Control | 5 | n/a |
| | <i>Or</i> | | |
| N/A* | Test sample | 5 | |
| | <i>Or</i> | | |
| N/A* | H ₂ O (molecular biology grade) | 5 | |

*5 µl of sample is recommended and should be used in most experiments. However, 2-8 µl of sample may be used. The volume of water must then be adjusted to ensure that the total reaction volume is 20 µl. Please add the same total amount of reagents plus DNA to all PCR vessels. The Positive Control should be used at 5 µl and the volume of additional water adjusted to ensure that the total reaction volume is 20 µl

**To compensate for any loss during pipetting it may be necessary to prepare an additional volume of reaction mix, a 5% excess is usually sufficient.

Note: Seal the PCR plate using PCR caps and centrifuge briefly. Every well should be sealed tightly to avoid evaporation. **4titude** 96-Well PCR Plate (cat. 4ti-0750/TA) and 8-Strip PCR Caps (cat. 4ti-0751) are recommended for good results.

7.4 Real-Time PCR instrument settings

The assay has been optimized for the BioRad CFX96 and the SLAN-96P Real-time PCR system. The sections of protocol describing run settings and data analysis parameters are specific for these systems.

7.5 Real-Time PCR instrument settings for Bio-Rad CFX96

- Place the plate in the instrument.
- PCR volume is set to 20 µl and “none” is selected for passive reference.
- Select **all channels** detection for all wells in use.
- The PCR run is performed using the Standard Run Mode with cycling conditions as described in the table below

| Stage | Cycles | Temperature (°C) | Duration | Data collection |
|------------------------------|--------|------------------|----------|---|
| Reverse Transcription | 1 | 55 | 10 min | |
| Hotstart | 1 | 95 | 3 min | |
| Amplification | 48 | 95 | 10 sec | |
| | | 60 | 30 sec | End-point point fluorescence collection |
| | | 69 | 15 sec | |
| Dissociation | 1 | 95 | 15 sec | |
| | | 25 | 30 sec | Real-time point fluorescence collection from 25°C to 75°C |
| | | 75 | 5 sec | |

7.5.1.1 Software for Data analysis

Perform data analysis using BioRad CFX Manager (Version 3.1).

7.5.1.2 Sample analysis

Amplification in the FAM, HEX and ROX channels indicates the presence of the viral E6/E7 mRNA targeted by the assay.

Amplification in Cy5 channel serves as an internal PCR reaction control for each sample. Every well should give a Cy5 signal (Ct<40) from transcripts of a human endogenous internal control (IC) gene. For NTC: No amplification (or Ct> 40) should be detected in FAM, HEX, ROX and CY5 channels.

7.5.1.3 Setting the baseline, threshold and Ct calling

The baseline should be set to a range that eliminates the background fluorescence found in the early cycles of amplification, but which does not overlap the area where amplifications signals rise above the background. Baseline may be set automatically if it gives a suitable value.

The cycle number at which a signal is detected above background fluorescence is termed the cycle threshold (Ct). Select the threshold for Ct determination as close as possible to the base of the exponential phase. As an indication, select a threshold that gives in all channels Ct values for PC between 25 and 33.

Obtain threshold cycle (Ct) values for each channel, the test is considered valid only if:

- NTC shows no amplification (or Ct >40) in FAM, HEX, ROX and CY5 channels.
- Positive Control shows amplification in FAM, HEX, ROX and Cy5 channels (Ct<35).
- Amplification of Internal Control in clinical sample is Ct<40 in Cy5 channel.

7.5.1.4 Viral identification via melting profile analysis


The term melting profile refers to the melting curves (per channel) generated during the dissociation stage of the reaction (from 25°C to 75°C). The melting profile obtained per channel is a combination of the melting data generated from all viral types included in this channel. Each viral type has a unique melting profile consisting of a unique melting temperature and shape. A change in this characteristic melting profile, in comparison with the NTC reference melting profile, shows the sample to be positive for the respective viral type(s).

The characteristic melting profile for each viral type are shown in the figures below:







Specimen genotyping results are interpreted as shown below.

- Y axis denotes derivative fluorescence
- X axis shows temperature
- Melting curve profiles of NTC that are suitable for analysis are shown below for each fluorescent channel
- Differences in melting curves that specifically identify viral types are indicated by arrows






First, check the melting curve profiles of the calibrator in Cy5 channel, the peak is usually detected at 43.0 ± 0.5 °C (see Figure below). If the peak has shifted to a lower or higher temperature, a similar shift will be observed in all the melting curves in all channels and so should be taken into consideration during the analysis.









Melting curves for viral types detected in FAM channel

| Genotyping result | Representative melting curve profile | Features |
|---|--|---------------|
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 16 (HPV16) |  | 48.0 ± 1.0 °C |
| Human Papillomavirus type 31 (HPV31) |  | 34.0 ± 1.0 °C |
| Human Papillomavirus type 66 (HPV66) |  | 55.0 ± 1.0 °C |
| Human Papillomavirus type 39 (HPV39) |  | 40.0 ± 1.0 °C |
| Human Papillomavirus type 59 (HPV59) |  | 59.0 ± 1.0 °C |

Melting curves for viral types detected in HEX channel

| Genotyping result | Representative melting curve profile | Features |
|---|--|--------------------------------|
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 56 (HPV56) |  | T _m : 58.0 ± 1.0 °C |
| Human Papillomavirus type 33 (HPV33) |  | T _m : 34.0 ± 1.0 °C |
| Human Papillomavirus type 45 (HPV45) |  | T _m : 53.0 ± 1.0 °C |
| Human Papillomavirus type 68 (HPV68) |  | T _m : 44.0 ± 1.0 °C |

Melting curves for viral types detected in ROX channel

| Genotyping result | Representative melting curve profile | Features |
|---|--|-------------------------------|
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 52 (HPV52) |  | T _m : 54.0 ± 1.0°C |
| Human Papillomavirus type 18 (HPV18) |  | T _m : 47.0 ± 1.0°C |
| Human Papillomavirus type 35 (HPV35) |  | T _m : 41.0 ± 1.0°C |
| Human Papillomavirus type 51 (HPV51) |  | T _m : 58.0 ± 1.0°C |
| Human Papillomavirus type 58 (HPV58) |  | T _m : 32.0 ± 1.0°C |

The table below gives information on sample results and suggested outcomes in different scenarios for the MPAHPV004, Papilloplex® HR-HPV mRNA kit.

| Amplification in Cy5 | +RT-PCR Amplification in FAM, HEX or ROX | Changes in melting profile | Sample result and suggested actions |
|----------------------|--|----------------------------|---|
| Yes | Yes | Yes | Presence of human and Virus mRNA |
| Yes | No | No | Viral E6/E7 mRNA Negative |
| Yes | Yes | No | Viral E6/E7 mRNA Positive; should be re-tested for genotyping results |
| Yes/No | No | Yes | Sample Invalid; should be re-tested |
| No | Yes | Yes | Sample Invalid; should be re-tested. If the same result is encountered again then the sample is considered positive; failure to amplify the internal control is due to partial degradation of the RNA in the sample |
| No | Yes | No | Sample Invalid; should be re-tested |
| No | No | No | Sample Invalid; should be re-tested |
| No | No | Yes | Sample Invalid; should be re-tested |

7.6 Real-Time PCR instrument settings for SLAN96P

- Place the plate in the instrument.
- PCR volume is set to 20 µl and “none” is selected for passive reference.
- Select **all channels** detection for all wells in use.
- The PCR run is performed using the **Melting Curve Mode** with cycling conditions as described in the table below

| Stage | Cycles | Temperature (°C) | Duration | Data collection |
|------------------------------|--------|------------------|----------|--|
| Reverse Transcription | 1 | 55 | 10 min | |
| Hotstart | 1 | 95 | 3 min | |
| Amplification | 48 | 95 | 10 sec | |
| | | 60 | 30 sec | End-point point fluorescence collection |
| | | 69 | 15 sec | |
| Dissociation | 1 | 95 | 15 sec | |
| | | 27 | 30 sec | Real-time point fluorescence collection every 0.5°C (stepwise) from 27°C to 65°C |

7.6.1.1 Software for data analysis

Perform data analysis using SLAN Real-time PCR System (Version 8.2.2).

7.6.1.2 Sample analysis

Amplification in the FAM, HEX and ROX channels indicates the presence of the viral DNA targeted by the assay.

Amplification in Cy5 channel serves as an internal PCR reaction control for each sample. Every tested samples should give a Cy5 signal ($Ct < 35$) from DNA of a human endogenous internal control (IC) gene. For NTC: No amplification (or $Ct > 38$) should be detected in FAM, HEX or ROX channel and no amplification or $Ct > 35$ should be detected in CY5 channel.

7.6.1.3 Setting the baseline, threshold and Ct calling

The baseline should be set to a range that eliminates the background fluorescence found in the early cycles of amplification, but which does not overlap the area where amplifications signals rise above the background. Baseline may be set automatically if it gives a suitable value.

The cycle number at which a signal is detected above background fluorescence is termed the cycle threshold (Ct). Select the threshold for Ct determination as close as possible to the base of the exponential phase. As an indication, select a threshold that gives in all channels Ct values for PC between 25 and 33.

Obtain threshold cycle (Ct) values for each channel, the test is considered valid only if:

- NTC shows no amplification (or $Ct > 40$) in FAM, HEX, ROX and CY5 channels.
- Positive Control shows amplification in FAM, HEX, ROX and Cy5 channels ($Ct < 35$).
- Amplification of Internal Control in clinical sample is $Ct \leq 40$ in Cy5 channel.

7.6.1.4 Viral identification via melting profile analysis


Before checking for the melting profile, **sample must show amplification with $Ct \leq 40$ to be considered as virus positive and considered for genotyping using the melt curve profile.** The term melting profile refers to the melting curves (per channel) generated during the dissociation stage of the reaction (from 25°C to 75°C). The melting profile obtained per channel is a combination of the melting data generated from all viral types included in this channel. Each viral type has a unique melting profile consisting of a unique melting temperature and shape. A change in this characteristic melting profile, in comparison with the NTC reference melting profile, shows the sample to be positive for the respective viral type(s).







Each viral type shows a characteristic melting profile which differs to the melting curve observed for the NTC in the given channel. The difference between the melting curve profiles indicates the viral type present in the sample as shown in the figures below.






Specimen genotyping results are interpreted as shown below.







- Y axis denotes derivative fluorescence
- X axis shows temperature
- Melting curve profiles of NTC that are suitable for analysis are shown below for each fluorescent channel
- Differences in melting curves that specifically identify viral types are shown by arrows

First, check the melting curve profiles of the calibrator in Cy5 channel, the peak is usually detected at 44.0 ± 0.5 °C. If the peak has shifted to a lower or higher temperature, a similar shift will be observed in all the melting curves in all channels and so should be taken into consideration during the analysis.



| Melting curves for viral types detected in FAM channel | | |
|--|--|---------------|
| Genotyping result | Representative melting curve profile | Features |
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 16 (HPV16) |  | 48.0 ± 1.0 °C |
| Human Papillomavirus type 31 (HPV31) |  | 35.0 ± 1.0 °C |
| Human Papillomavirus type 66 (HPV66) |  | 55.0 ± 1.0 °C |
| Human Papillomavirus type 39 (HPV39) |  | 40.0 ± 1.0 °C |
| Human Papillomavirus type 59 (HPV59) |  | 60.0 ± 1.0 °C |

| Melting curves for viral types detected in HEX channel | | |
|--|--|-------------------|
| Genotyping result | Representative melting curve profile | Features |
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 56 (HPV56) |  | Tm: 59.0 ± 1.0 °C |
| Human Papillomavirus type 33 (HPV33) |  | Tm: 34.0 ± 1.0°C |
| Human Papillomavirus type 45 (HPV45) |  | Tm: 53.0 ± 1.0°C |
| Human Papillomavirus type 68 (HPV68) |  | Tm: 45.0 ± 1.0°C |

| Melting curves for viral types detected in ROX channel | | |
|--|--|------------------|
| Genotyping result | Representative melting curve profile | Features |
| Standard/expected NTC melting curve profile |  | N/A |
| Human Papillomavirus type 52 (HPV52) |  | Tm: 54.0 ± 1.0°C |
| Human Papillomavirus type 18 (HPV18) |  | Tm: 48.0 ± 1.0°C |
| Human Papillomavirus type 35 (HPV35) |  | Tm: 42.0 ± 1.0°C |
| Human Papillomavirus type 51 (HPV51) |  | Tm: 59.0 ± 1.0°C |
| Human Papillomavirus type 58 (HPV58) |  | Tm: 34.0 ± 1.0°C |

The table below gives information on sample results and suggested outcomes in different scenarios for the MPAHPV004 Papilloplex® HR-HPV mRNA kit.

| Amplification in Cy5 (Ct ≤ 40) | +RT-PCR Amplification in FAM, HEX or ROX (Ct ≤ 40) | Changes in melting profile | Sample result and suggested actions |
|--------------------------------|--|----------------------------|---|
| Yes | Yes | Yes | Presence of human and Virus mRNA |
| Yes | No | No | Viral E6/E7 mRNA Negative |
| Yes | Yes | No | Viral E6/E7 mRNA Positive; should be re-tested for genotyping results |
| Yes/No | No | Yes | Sample Invalid; should be re-tested |
| No | Yes | Yes | Sample Invalid; should be re-tested. If the same result is encountered again then the sample is considered positive; failure to amplify the internal control is due to partial degradation of the RNA in the sample |
| No | Yes | No | Sample Invalid; should be re-tested |
| No | No | No | Sample Invalid; should be re-tested |
| No | No | Yes | Sample Invalid; should be re-tested |

8 Troubleshooting

Should you encounter problems please consult the table below:

| Observation | Probable Cause | Solution |
|--|--|--|
| Absence of amplification in test samples | Presence of PCR inhibitors | The performance of the kit may be adversely affected by known PCR inhibitors co-extracted from patient samples (such as blood, excessive mucus, pharmaceutical preparations such as lubricant gels, spermicide creams etc.) We suggest repeating the test or obtaining a new patient sample/nucleic acid extraction. |
| | Insufficient nucleic acid in test sample | Repeat processing of the same sample using a greater test volume or obtain a new patient sample. |
| Absence of amplification in positive control | Instrument faulty | Check the instrument calibration records and confirm it is working. |
| | Kit stored at wrong temperature or under wrong conditions | Check the storage temperature and whether the contents were exposed to prolonged direct sunlight. Also ensure the reagents are kept on ice during use and avoid excessive vortexing. |
| | Incorrect PCR cycling parameters | Verify that PCR cycling parameters correspond to those recommended above. |
| Melting profile that does not resemble the reference NTC melt profile | Two adjacent changes in melting curves can make it appear different from the reference profiles | Two or more viral types might be present in the sample. |
| | Unequal volume of liquid in tubes due to evaporation of liquid from one tube or pipetting errors | Proceed with caution. If the melting curves are not possible to visually align, repeat the samples and NTC. |

9 Performance

9.1 Stability

The purpose of stability testing is to provide evidence on how the quality of Papilloplex® HR-HPV mRNA kit varies under the different environmental factors such as temperature or freeze-thawing.

could undergo 4 cycles of thaw-freeze without any loss in performance.

Papilloplex® HR-HPV mRNA kit could go through 8h at 37°C, 12h at room temperature and freeze again at -20°C without any loss in performance.

Papilloplex® HR-HPV mRNA kit could be set-up at room temperature up to 30°C without any loss in performance.

After setting-up the Papilloplex® HR-HPV mRNA kit, a delay in starting the PCR reaction (2 hours 40 min) does not affect the performance of the kit.

9.2 Limit of Detection (LOD)

The Limit of detection of Papilloplex® HR-HPV mRNA kit was determined on the Bio-Rad CFX96 real-time PCR instrument using quantified DNA plasmids containing E6/E7 partial sequences for each HR-HPV targeted by the assay.

| Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|--------|--|--------------|-----------------------|------|
| HPV16 | 200 000 | 24.2 | 0.4 | 1.5 |
| | 10 000 | 27.1 | 0.4 | 1.4 |
| | 500 | 31.5 | 0.3 | 1.0 |
| | 100 | 34.2 | 0.5 | 1.4 |
| HPV18 | 200 000 | 23.5 | 0.5 | 2.2 |
| | 10 000 | 25.5 | 0.3 | 1.4 |
| | 500 | 30.3 | 0.3 | 1.1 |
| | 100 | 32.9 | 0.3 | 1.0 |
| HPV31 | 200 000 | 23.2 | 0.4 | 1.9 |
| | 10 000 | 27.2 | 0.8 | 2.9 |
| | 500 | 31.4 | 0.6 | 1.8 |
| | 100 | 34.3 | 0.8 | 2.3 |
| HPV33 | 200 000 | 23.8 | 0.3 | 1.1 |
| | 10 000 | 27.7 | 0.4 | 1.6 |
| | 500 | 31.9 | 0.6 | 1.9 |
| | 100 | 34.7 | 0.7 | 2.1 |
| HPV35 | 200 000 | 23.5 | 0.4 | 1.9 |
| | 10 000 | 27.2 | 0.6 | 2.0 |
| | 500 | 31.9 | 0.5 | 1.4 |
| | 100 | 34.6 | 0.3 | 0.9 |
| HPV39 | 200 000 | 21.9 | 0.5 | 2.4 |
| | 10 000 | 27.0 | 0.5 | 1.8 |
| | 500 | 31.7 | 0.5 | 1.6 |
| | 100 | 34.7 | 0.8 | 2.4 |
| HPV45 | 200 000 | 20.4 | 0.1 | 0.6 |

| Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|--------------|--|--------------|-----------------------|------|
| | 10 000 | 26.1 | 0.5 | 1.9 |
| | 500 | 30.3 | 0.3 | 1.0 |
| | 100 | 32.7 | 0.5 | 1.5 |
| HPV51 | 200 000 | 23.4 | 0.5 | 2.0 |
| | 10 000 | 29.8 | 0.6 | 0.6 |
| | 500 | 34.2 | 0.6 | 1.9 |
| | 100 | 37.6 | 0.7 | 1.9 |
| HPV52 | 200 000 | 23.8 | 0.4 | 1.9 |
| | 10 000 | 26.7 | 0.5 | 1.9 |
| | 500 | 31.7 | 0.3 | 0.8 |
| | 100 | 34.4 | 0.5 | 1.4 |
| HPV56 | 200 000 | 25.2 | 0.9 | 3.5 |
| | 10 000 | 27.1 | 0.7 | 2.4 |
| | 500 | 32.3 | 0.3 | 1.1 |
| | 100 | 35.1 | 0.5 | 1.5 |
| HPV58 | 200 000 | 19.7 | 0.3 | 1.3 |
| | 10 000 | 25.4 | 0.4 | 1.7 |
| | 500 | 29.5 | 0.4 | 1.3 |
| | 100 | 31.9 | 0.4 | 1.4 |
| HPV59 | 200 000 | 23.3 | 0.3 | 1.4 |
| | 10 000 | 28.2 | 0.5 | 1.9 |
| | 500 | 33.5 | 1.5 | 4.4 |
| | 100 | 36.5 | 0.8 | 2.1 |
| HPV66 | 200 000 | 26.7 | 0.3 | 1.1 |
| | 10 000 | 29.5 | 0.7 | 2.3 |
| | 500 | 34.8 | 0.6 | 1.7 |
| | 100 | 38.0 | 0.7 | 1.7 |
| HPV68 | 200 000 | 22.1 | 0.3 | 1.1 |
| | 10 000 | 26.5 | 0.8 | 2.9 |
| | 500 | 31.3 | 0.5 | 1.6 |
| | 100 | 34.3 | 0.5 | 1.6 |

9.3 Analytical Specificity (cross-reactivity with others HPV types)

The BLAST (<https://blast.ncbi.nlm.nih.gov/Blast.cgi>).search showed that the primers and the probes included in the Papilloplex® HR-HPV mRNA kit do not cross-react with all sequences known for low-risk (LR) HPV types. Hence all sequences of known LR-HPV types are included in NCBI databases, finding no hit LR-HPV sequence with the potential of producing a detectable PCR-fragment presents a strong evidence that none of the LR-HPV types will give a false-positive signal.

9.4 Analytical Specificity (cross-reactivity with other microorganisms)

Potential cross-reactivity with non-HPV microorganisms reasonably expected to be present at the site of clinical sample collection was assessed, using quantified genomic DNA samples at 10 000 copies per reaction of *Lactobacillus acidophilus*, *Candida albicans*, *Mycoplasma hominis*, *Trichomonas vaginalis C-1*, *Garnerella vaginalis*, *Staphylococcus aureus*, *Chlamydia trachomatis Serovar E*, *Human herpes virus 2*, *Escherichia coli*, and *Neisseria gonorrhoeae*. No cross-reactivity were observed.

9.5 Precision

9.6 Precision Intra-Assay

The precision Intra-assay was performed by one operator on one instrument using one batch of Papilloplex® HR-HPV mRNA kit. The study consisted of 3 runs per instrument and 3 concentrations were assessed : Low (10 times LOD; 1000 copies/reaction), Medium(100 times the LOD; 10 000 copies/reaction) and high (>100 times the LOD; 100 000 copies per reaction). The results are summarised below:

| Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|--------|--|--------------|-----------------------|------|
| HPV16 | 100 000 | 25.2 | 0.4 | 1.8 |
| | 10 000 | 27.1 | 0.3 | 1.3 |
| | 1 000 | 30.9 | 0.4 | 1.2 |
| HPV18 | 100 000 | 23.7 | 0.3 | 1.4 |
| | 10 000 | 27.1 | 0.4 | 1.7 |
| | 1 000 | 29.8 | 0.4 | 1.4 |
| HPV31 | 100 000 | 25.9 | 1.3 | 4.9 |
| | 10 000 | 26.0 | 0.4 | 1.4 |
| | 1 000 | 32.9 | 0.4 | 1.1 |
| HPV33 | 100 000 | 26.5 | 0.3 | 1.2 |
| | 10 000 | 26.8 | 0.4 | 1.6 |
| | 1 000 | 33.7 | 0.3 | 0.8 |
| HPV35 | 100 000 | 24.9 | 0.2 | 0.8 |
| | 10 000 | 27.7 | 0.4 | 1.3 |
| | 1 000 | 32.6 | 0.4 | 1.3 |
| HPV39 | 100 000 | 25.0 | 0.3 | 1.3 |
| | 10 000 | 28.7 | 0.3 | 1.2 |
| | 1 000 | 29.6 | 0.8 | 2.8 |
| HPV45 | 100 000 | 24.5 | 0.5 | 1.9 |
| | 10 000 | 25.3 | 0.2 | 0.9 |
| | 1 000 | 28.7 | 0.3 | 1.2 |
| HPV51 | 100 000 | 27.4 | 0.3 | 1.2 |
| | 10 000 | 31.3 | 0.2 | 0.6 |
| | 1 000 | 32.1 | 0.9 | 2.9 |
| HPV52 | 100 000 | 25.0 | 0.4 | 1.7 |
| | 10 000 | 26.8 | 0.3 | 1.3 |
| | 1 000 | 30.9 | 0.4 | 1.3 |
| HPV56 | 100 000 | 25.3 | 0.3 | 1.2 |
| | 10 000 | 28.8 | 0.4 | 1.2 |
| | 1 000 | 31.7 | 0.4 | 1.3 |
| HPV58 | 100 000 | 24.4 | 0.4 | 1.6 |
| | 10 000 | 24.7 | 0.3 | 1.2 |
| | 1 000 | 28.4 | 0.3 | 1.0 |
| HPV59 | 100 000 | 25.3 | 0.4 | 1.7 |
| | 10 000 | 29.4 | 0.2 | 0.6 |
| | 1 000 | 31.5 | 0.5 | 1.6 |

| Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|--------|--|--------------|-----------------------|------|
| HPV66 | 100 000 | 27.6 | 0.6 | 2.0 |
| | 10 000 | 30.1 | 0.4 | 1.2 |
| | 1 000 | 35.1 | 0.6 | 1.7 |
| HPV68 | 100 000 | 24.0 | 0.2 | 1.0 |
| | 10 000 | 27.5 | 0.4 | 1.5 |
| | 1 000 | 29.7 | 0.4 | 1.4 |

9.7 Precision Inter-Assay

Verification of Reproducibility was performed over 8-days, three time per day by two operators on two instruments using three batches of Papilloplex® HR-HPV mRNA kit. The study consisted of 24 runs per platform (Bio-rad CFX96 or Slan96P) and 3 concentrations were assessed : Low (10 times LOD; 1000 copies/reaction), Medium(100 times the LOD; 10 000 copies/reaction) and high (>100 times the LOD; 100 000 copies per reaction). No significant variability between operators, runs and batches was recorded at any concentration tested.

The results are summarised below:

| Instruments | Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|--------------|---------|--|--------------|-----------------------|------|
| BioRad CFX96 | HPV16 | 100 000 | 24.8 | 1.4 | 5.6 |
| | | 10 000 | 27.1 | 0.7 | 2.4 |
| | | 1 000 | 30.6 | 1.1 | 3.4 |
| | HPV18 | 100 000 | 23.0 | 0.8 | 3.3 |
| | | 10 000 | 26.3 | 1.0 | 3.9 |
| | | 1 000 | 30.3 | 1.3 | 4.1 |
| | HPV31 | 100 000 | 23.3 | 1.0 | 4.3 |
| | | 10 000 | 27.0 | 1.4 | 5.3 |
| | | 1 000 | 31.1 | 1.7 | 5.6 |
| | HPV33 | 100 000 | 23.5 | 1.0 | 4.4 |
| | | 10 000 | 27.1 | 1.5 | 5.5 |
| | | 1 000 | 31.2 | 1.5 | 4.9 |
| | HPV35 | 100 000 | 24.3 | 0.8 | 3.1 |
| | | 10 000 | 28.0 | 0.7 | 2.4 |
| | | 1 000 | 30.8 | 0.8 | 2.5 |
| | HPV39 | 100 000 | 23.3 | 0.7 | 3.2 |
| | | 10 000 | 26.9 | 1.2 | 4.5 |
| | | 1 000 | 30.1 | 0.9 | 2.9 |
| | HPV45 | 100 000 | 22.2 | 1.4 | 6.1 |
| | | 10 000 | 25.1 | 0.8 | 3.0 |
| | | 1 000 | 29.3 | 1.6 | 5.5 |
| | HPV51 | 100 000 | 25.6 | 1.0 | 4.0 |
| | | 10 000 | 29.7 | 1.2 | 4.0 |
| | | 1 000 | 33.1 | 1.1 | 3.4 |
| HPV52 | 100 000 | 24.3 | 1.6 | 6.4 | |

| Instruments | Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|-------------|--------|--|--------------|-----------------------|------|
| Slan 96P | | 10 000 | 27.0 | 0.8 | 2.8 |
| | | 1 000 | 30.9 | 0.9 | 2.9 |
| | HPV56 | 100 000 | 25.1 | 1.1 | 4.4 |
| | | 10 000 | 28.5 | 1.5 | 5.2 |
| | | 1 000 | 32.8 | 1.6 | 5.0 |
| | | HPV58 | 100 000 | 22.1 | 0.8 |
| | 10 000 | | 25.2 | 0.8 | 1.6 |
| | | 1 000 | 29.3 | 1.6 | 5.3 |
| | | HPV59 | 100 000 | 23.8 | 0.8 |
| | 10 000 | | 27.5 | 1.3 | 4.6 |
| | | 1 000 | 32.0 | 1.1 | 3.3 |
| | | HPV66 | 100 000 | 27.1 | 0.8 |
| | 10 000 | | 30.9 | 0.8 | 2.7 |
| | | 1 000 | 33.7 | 1.0 | 2.9 |
| | | HPV68 | 100 000 | 22.6 | 0.7 |
| | 10 000 | | 26.2 | 1.1 | 4.3 |
| | | 1 000 | 30.2 | 0.9 | 3.1 |
| | | HPV16 | 100 000 | 21.8 | 0.5 |
| | 10 000 | | 26.2 | 1.2 | 4.8 |
| | | 1 000 | 29.3 | 1.2 | 4.2 |
| | | HPV18 | 100 000 | 18.6 | 0.9 |
| | 10 000 | | 22.9 | 1.2 | 5.2 |
| | | 1 000 | 26.0 | 1.1 | 4.2 |
| | | HPV31 | 100 000 | 21.2 | 0.5 |
| | 10 000 | | 26.3 | 1.8 | 6.8 |
| | | 1 000 | 28.8 | 1.6 | 5.7 |
| | | HPV33 | 100 000 | 21.0 | 0.5 |
| | 10 000 | | 25.9 | 1.8 | 7.1 |
| | | 1 000 | 28.4 | 1.7 | 6.0 |
| | | HPV35 | 100 000 | 20.8 | 0.5 |
| | 10 000 | | 25.1 | 1.0 | 3.8 |
| | | 1 000 | 28.5 | 0.7 | 2.4 |
| | | HPV39 | 100 000 | 21.5 | 0.8 |
| | 10 000 | | 26.1 | 1.5 | 5.7 |
| | | 1 000 | 29.3 | 1.7 | 5.7 |
| | | HPV45 | 100 000 | 19.3 | 0.5 |
| 10 000 | 24.0 | | 1.3 | 5.6 | |
| | 1 000 | 27.1 | 1.3 | 4.7 | |
| | HPV51 | 100 000 | 20.4 | 0.7 | 3.5 |
| 10 000 | | 25.5 | 1.6 | 6.1 | |
| | 1 000 | 29.1 | 1.7 | 5.7 | |
| | HPV52 | 100 000 | 19.5 | 0.6 | 3.2 |
| 10 000 | | 24.4 | 1.6 | 6.5 | |

| Instruments | Target | Concentration Copies per reactions | Ct (Mean) | Standard Deviation | CV % |
|-------------|--------|--|--------------|-----------------------|------|
| | | 1 000 | 27.9 | 1.0 | 3.7 |
| | HPV56 | 100 000 | 22.4 | 0.8 | 3.7 |
| | | 10 000 | 27.0 | 1.7 | 6.4 |
| | | 1 000 | 30.2 | 1.5 | 4.8 |
| | HPV58 | 100 000 | 18.3 | 0.5 | 3.0 |
| | | 10 000 | 23.2 | 1.2 | 5.2 |
| | | 1 000 | 26.2 | 1.3 | 5.0 |
| | HPV59 | 100 000 | 22.6 | 0.8 | 3.5 |
| | | 10 000 | 28.3 | 1.5 | 5.4 |
| | | 1 000 | 30.5 | 1.1 | 3.6 |
| | HPV66 | 100 000 | 27.7 | 1.0 | 3.7 |
| | | 10 000 | 32.2 | 2.3 | 7.1 |
| | | 1 000 | 35.4 | 2.2 | 6.3 |
| | HPV68 | 100 000 | 19.7 | 0.4 | 1.8 |
| | | 10 000 | 25.1 | 1.3 | 5.2 |
| | | 1 000 | 27.2 | 0.7 | 2.7 |

9.8 Clinical Performance

9.9 Clinical specimens

225 samples (LBC in PreserCyt medium (ThinPrep®) were obtained from the Scottish HPV reference laboratory. Nucleic acids were extracted using either Quick DNA/RNA Viral Kit from Zymo Research or Viral-PrepAdem-Kit from Ademtech.

A summary of the distribution of the samples according to their Cytology and Histology diagnoses is presented below:

| Cytology/Histology | Not available | Cytology Neg so no Histology | Negative for dysplasia | Koilocytosis | CIN1 | CIN2 | CIN3 |
|--------------------------|---------------|------------------------------|------------------------|--------------|------|------|------|
| Negative | 0 | 85 | 0 | 0 | 0 | 0 | 0 |
| Bordeline Squamous | 29 | 0 | 6 | 1 | 3 | 1 | 0 |
| Low grade dysk | 3 | 0 | 10 | 1 | 2 | 0 | 4 |
| High Grade moderate dysk | 0 | 0 | 6 | 1 | 3 | 6 | 5 |
| High Grade severe dysk | 0 | 0 | 1 | 0 | 2 | 15 | 41 |

9.10 Positive rate of HPV E6/E7 mRNA assay by histological diagnosis

| | Negative for dysplasia | Koilocytosis | CIN1 | CIN2 | CIN3 |
|---------------|------------------------|--------------|------|------|------|
| m RNA+ | 14 | 1 | 9 | 19 | 47 |
| mRNA - | 7 | 2 | 1 | 3 | 3 |
| Positive rate | 67% | 33% | 90% | 86% | 94% |

9.11 Concordance between Papilloplex® HR-HPV mRNA kit and Cytology diagnosis

| | GeneFirst | Cytology | GeneFirst | Cytology |
|------------------------------------|-----------|------------|-----------|--------------|
| Histology | E6/E7- | <Low Grade | E6/E7+ | > High Grade |
| < CIN1 | 10 | 21 | 24 | 13 |
| > CIN2 + | 6 | 5 | 66 | 67 |
| > CIN3+ | 3 | 4 | 47 | 46 |
| concordance in CIN 2+ (86%) | | | | |
| concordance in CIN3+ (85%) | | | | |

9.12 Correlation of Papilloplex® HR-HPV mRNA kit or Cytology with Histological diagnoses (CIN2+)

| | Sensitivity % (CI) | Specificity % (CI) | Predictive Positive Value (PPV) % (CI) | Negative Predictive Value (NPV) % (CI) | Accuracy % (CI) |
|------------------|-------------------------|-------------------------|---|---|------------------------|
| Cytology | 93.1 % (84.5 - 97.7) | 61.8 % (43.6 - 77.8) | 83.8 % (77.0 - 88.8) | 80.8 % (63.4 - 91.1) | 83.0% (74.5- 89.6) |
| GeneFirst | 91.7% (82.7 -96.9) | 70.6 % (52.5 - 84.9) | 80% (64.3 -89.9) | 86.8 % (79.6 -91.8) | 84.9% (76.7 - 91.1) |

10 Product Specification

| | |
|--|---|
| Technology | Reverse-Transcription Real-Time PCR |
| Target Sequence | E6/E7 region |
| Clinical Specificity | 70.6 % |
| Clinical Sensitivity | 91.7% |
| Sensitivity (LOD) | 100 copies per reactions |
| Kit Storage | -20 °C |
| Sample Material | No Sample material supplied with the Kit, however Genefirst recommends the following DNA Extraction methods: <ul style="list-style-type: none"> - Quick DNA/RNA Viral Kit from Zymo Research - Viral-PrepAdem-Kit from Ademtech |
| Validated Real-Time PCR Devices | Biorad-CFX96 and SLAN96P |
| Quality Control | Positive and internal cellular controls of PCR amplification and sample integrity |
| Certification | CE IVD for <i>in vitro</i> Diagnostics Use |

11 Symbols

| | | | |
|---|---|---|------------------------------------|
|  | Consult instructions for use |  | Upper limit of temperature -20°C |
|  | Catalogue number |  | Keep away from sunlight |
|  | Date of manufacture |  | Batch code |
|  | Manufacturer |  | In vitro diagnostic medical device |
|  | Use-by-date |  | Do not use if package damaged |
|  | Authorised representative in the European Community | | |

Please Note: Information in these instructions for use is subject to change without notice and does not represent commitment on the part of GeneFirst. No part of these instructions for use may be reproduced or transmitted in any form or by any means, electronic or mechanical, including photocopying and recording for any purpose without the written permission of GeneFirst.

12 Customer contact information

For all sales order processing, training and technical support enquiries, please contact the following:



GeneFirst Limited

Unit 2 The Quadrant,
Abingdon Science Park,
Abingdon,
Oxfordshire.
OX14 3YS
UK

Customer Service & Sales Enquiries:

Telephone: +44 (0)1865 407 400

Email: sales@genefirst.com

www.genefirst.com



Advena Ltd. Tower Business Centre, 2nd Flr.,
Tower Street, Swatar, BKR 4013 Malta

Part no. IFU0709

Issue 4.0

March 2022

© GeneFirst 2022