# Istruzioni d'uso

ver. 8 - 26/05/2023

Biologia Molecolare

## CVD14

**REF**

AC084

20 TEST

CND

W0106010499

**CE**

**IVD**

**NUCLEAR LASER MEDICINE S.r.l.**

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
**IVD**

**NUCLEAR LASER MEDICINE S.r.l.**

HEADQUARTER: Via Cascina Conighetto – BUSINESS OFFICES: Viale delle Industrie, 3 - 20049 SETTALA (MI), ITALY

Phone: (+39) 02/95.24.51 - Fax (+39) 02/95.24.52.37 - (+39) 02/95.24.52.38 WEB: [www.nlm.it](http://www.nlm.it) - E-MAIL: [segreteria@nlm.it](mailto:segreteria@nlm.it)

*ISO 9001 Quality Management System Organization certified and ISO 13485 Medical Sector Quality Management System certified (IMQ certification body - CSQ certification).*

## INTENDED USE

The **CVD 14** device provides reagents for the identification of the main gene mutations and polymorphisms involved in arterial and deep venous thrombosis (Factor V R506Q and H1299R, Factor II G20210A, MTHFR C677T and A1298C, CBS 844ins68, PAI-1, ACE, AGT, GPIIIa HPA-1 a/b, ATR-1,  $\beta$ - fibrinogen and Factor XIII) and in the cholesterol metabolism (ApoE), through the amplification of target sequences, reverse-hybridization and color development.

Amplification enzyme assigned: **orange** cap REF BA092, followed by size indication, which could change depending on enzyme quantity supplied (i.e.: BA092/50, BA092/60, etc).

This device has been validated with:

- a) manual column-based extraction (REF AA1001) starting from fresh blood or dried blood-spot and with the most common automatic systems for nucleic acids purification.
- b) PTC 100 (MJ Research); 2720 (Applied Biosystems), Mastercycler (Eppendorf); T100, C1000 and Mycycler (BioRad) thermal cycler instruments.
- c) automatic detection instruments: Tecan ProfiBlot 48 and Dynex Dynablot Heat

## INTRODUCTION

The development of cardiovascular disease in humans is the most common cause of morbidity and mortality in the Western world. CVD can be classified in broad terms as flow and pressure disorders of the cardiovascular system. Established risk factors include age, gender, diabetes mellitus, obesity, high serum cholesterol levels, hypertension, smoking and physical inactivity. These factors, however, do not explain all CVD cases. Infact in some cases the only evident risk factor is a positive familiar history of premature cardiovascular disease, implying genetic predisposition to the development of the disease. Therefore, in the last years a rapidly growing interest has been developed towards the genetic markers associated with cardiovascular disease risk factors, in order to find new prevention and/or therapeutic strategies. Ongoing studies are focusing on the possible interactions between the environmental and the genetic risk factors. Among the environmental components there are life-style (alcohol, smoking and drugs consumption), diet and stress; genetic susceptibility can be due to mutation and/or polymorphisms in several genes involved in the coagulation cascade, in the blood pressure regulation and in the metabolism of the lipids, glucose, homocysteine and iron. Among the candidate markers for inherited susceptibility to CVD the following genes have been selected: blood coagulation factors V (FV), II (prothrombin) and XIII (FXIII), methylenetetrahydrofolate reductase (MTHFR), cystathionine beta synthase (CBS), plasminogen activator inhibitor-1 (PAI-1), platelet glycoprotein IIIa (GPIIIa), angiotensin-converting enzyme (ACE), apolipoprotein E (Apo E), angiotensinogen (AGT), angiotensin II type 1 receptor (ATR-1) and  $\beta$ -Fibrinogen (FGB).

## REFERENCES

- Faioni E.M. et al. Blood 1999; 94 (9): 3062-6.
- Doggen et al. Circulation 1998 (97): 1037-1041.
- De Stefano et al. Blood 1998; 91 (10): 3562-3565.
- Wald et al. BMJ 2002 (325): 1202-1208.
- Gemmati et al. Haematologica 1999 (84): 824-828.
- Eriksson et al. Art Thromb and Vasc Biology 1998 (18): 20-26.
- Eichner et al. Am. J. Epidemiol. 2002 (155): 487-95.
- Weiss et al. N Engl J Med 1996 (334): 1090-4.
- Morise et al. J Intern Med 1995; 237 (2): 175-80.
- Crisan et al. J Mol Diagn 2000 (2): 105-115.
- Buraczyńska M et al. Kardiol Pol. 2003 Jan;58(1):1-9
- Blake G.J. et al. Eur Heart J 2001; 22 (24): 2262-2266.
- Gemmati et al. Wound Rep Reg 2004; 12 (5): 512-517.

## PRINCIPLE OF THE PROCEDURE

The test is based on the reverse-hybridization principle, where specific oligonucleotide probes immobilized as parallel lines on membrane-based strips hybridize with biotinylated PCR products. The exact match between probes and amplified product generates a signal exploiting the bond between biotin and streptavidin conjugated with alkaline phosphatase and a subsequent color developer.

The analysis consists of three steps:

1. DNA extraction from non-coagulated blood (see the procedure in the following *Extraction* paragraph).
2. PCR of target DNA sequences using biotinylated primers.
3. Hybridization of the amplification products with the oligonucleotide probes and colorimetric detection.

The device allows simultaneous detection of the following selected mutations/ polymorphisms:

1. Fattore V Leiden G1691A (R506Q)
2. Fattore V H1299R (HR2 haplotype)
3. Fattore II G20210A (Prothrombin)
4. MTHFR C677T
5. MTHFR A1298C
6. CBS 844ins68 I/D
7. PAI-1 (4G/5G)
8. GPIIIa T1565C (HPA-1 a/b)
9. ACE I/D
10. ApoE T112C\*
11. ApoE T158C\*
12. AGT g.9543T>C (M235T)
13. ATR-1 A1166C
14. FGB -455 G>A
15. Fattore XIII g.7130G>T (V34L)

⇒ \*ApoE: the following table shows the correspondence between the official nomenclature and the description of the investigated variants used in the kit:

Polymorphism	Nucleotide	Amino Acid	Description in the kit
rs429358	334T>C	C112R	T112C
rs7412	472C>T	R158C	T158C

## PRODUCT COMPOSITION

### BOX A (store at -25/-15°C)

REAGENTS	REF	Quantity	N° Vial
<b>PCR Mix 1</b> Solution containing oligonucleotides 1	KA832-n	1300 µl	1
<b>PCR Mix 2</b> Solution containing oligonucleotides 2	KA833-n	1300 µl	1

**BOX B (store at +2/+8°C)**

Reagents		REF	Quantity	N° Vial	
<b>Strips</b> Nitrocellulose membrane coated with oligonucleotides		-	KC021-n	20	1
<b>DNAT</b> Denaturation solution containing NaOH	 DANGER	H314, P280, P301+P330+P331, P303+P361+P353, P305+P351+P338	KA107-n	1,6 ml	1
<b>Hybridization (*)</b> Saline Solution containing preservatives	 WARNING	H319, EUH208, P280, P305+P351+P338, P337+P313	KA102-n	48 ml	1
<b>Wash Solution A</b> Saline solution containing detergent ad preservatives	-	-	KA103-n	136 ml	1
<b>Conjugate</b> Solution containing streptavidin labeled with alkaline phosphatase, stabilizers and preservatives	-	-	KA574-n	48 ml	1
<b>Wash Solution B</b> Saline Solution containing preservatives	-	EUH208	KA105-n	136 ml	1
<b>Color developer</b> Solution containing 5-Bromo-4-chloro-3-indolyl-phosphate and 4-Nitroblue tetrazolium	-	-	AA564-n	48 ml	1
<b>Transparent interpretation card</b>	-	-		1	
<b>Collector sheet</b> For strip storage, not suitable for results interpretation and/or as medical report	-	-		1	

(\*) The exemption from the requirements for labeling and packaging is determined by quantity and severity of danger reagent. The criteria are: **amounts <125 ml and categories not serious (ref. Regulation (EC) N ° 1272/2008 - list 1.5.2.1.1).**

The components volume show above refers to the standard size of the kit. Reduced packages of the kit are available on request for evaluation and / or demonstration of the product.

**STABILITY AND STORAGE**

- All opened or unopened reagents are stable up to the expiry date indicated on the label when stored at the correct temperature (**BOX A -25/-15°C; BOX B +2/+8°C**).
- At the end of each assay store reagents at the appropriate temperature.
- The kit should be kept isolated from any source of contaminating DNA, especially amplified products.
- Developed dry strips should be stored in the dark at room temperature (+15/+30°C).
- The vial containing the DNAT should be closed immediately after use.

## PRECAUTIONS

- Only professional and opportunely trained personnel should use this kit handle this product according to established good laboratory practices and universal precautions; wear personal protective apparel.
- When handling DNA samples or amplification products, use aerosol-resistant pipette tips to avoid pipettes contamination.
- Discard used materials as bio hazardous waste.
- Plastic trays can be reused maximum twice, if washed with water at the end of detection.
- Do not eat, drink, smoke or apply cosmetics in areas where reagents or specimens are handled.
- If skin or mucous membrane exposure occurs, immediately wash the area with copious amount of water. Seek medical advice.
- Do not use components beyond the expiration date.
- Test strips stored at +2/+8°C in the dark are stable until the expiry date.
- Do not mix reagents from different lots.
- Reagents have to be preserved separated from possible contaminants (as DNA samples and amplification products).
- It is recommended to work in three separated areas:
  - Area 1: pre-PCR (samples handling and extraction).
  - Area 2: Master Mix preparation.
  - Area 3: post-PCR (amplification and detection).
- Don't use the device if the box is damaged; contact the supplier.
- **It is advisable to have constant and uniform laboratory temperature, avoid to place the instruments near heating/cooling sources that may compromise the correct working.**
- **Use the transparent interpretation card of the specific lot.**
- **Hazard statements**
  - **H314:** Causes severe skin burns and eye damage.
  - **H319:** Causes serious eye irritation.
- **Precautionary statements**
  - **P280:** Wear protective gloves/protective clothing/eye protection/face protection.
  - **P301+P330+P331:** IF SWALLOWED: rinse mouth. Do NOT induce vomiting.
  - **P303+P361+P353:** IF ON SKIN (or hair): Remove/Take off immediately all contaminated clothing. Rinse skin with water/shower.
  - **P305+P351+P338:** IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
  - **P337+P313:** If eye irritation persists: Get medical advice/attention.
- **Supplemental hazard informations**
  - **EUH208:** It Contains CMIT/MIT. May cause allergic reaction (skin).

## MATERIALS REQUIRED BUT NOT SUPPLIED

### AREA 1

- DNA Extraction Kit
- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips and 1,5 ml tubes

### AREA 2

- Vertical downflow air box
- Dedicated adjustable volume pipettes set
- Aerosol-resistant pipettes tips
- 0,2 ml PCR tubes
- DNA polymerase with orange cap

### AREA 3

- Validated thermal cycler or equivalent instrument\*
- Agarose gel electrophoresis equipment (optional)
- Dedicated adjustable volume pipettes set and tips
- Clean tweezers to handle strips
- Plastic supports (trays) for manual or automatic detection

#### Manual strip processing:

- Water bath with inclined lid, shaking platform (50-100 rpm) and adjustable temperature ( $45 \pm 0,5$  °C)
- Thermometer
- Vacuum aspiration apparatus
- Deionized or distilled water
- Timer

#### Automatic strip processing:

- Tecan ProfiBlot 48, Dynex Dynablot Heat or equivalent instrument.

\* **WARNING:** if thermal cyclers different from those mentioned above are used, verify that the ramp rate (heating/cooling) is  $\leq 3^\circ\text{C}/\text{sec}$ . For T100 and C1000 (BioRad) set ramp rate as  $2,5^\circ\text{C}/\text{sec}$ .

## PROCEDURE

### DNA EXTRACTION

Manual purification: REF AA1001.

For extraction with automatic systems follow the specific manufacturer's IFUs.

Any other type of DNA extraction with similar results of concentration and purification can be used.

Use EDTA or citrate only as anticoagulant, do not use heparin. It is possible to use fresh blood, blood stored at +2/+8°C for no more than 2 days or stored at -25/-15°C for longer storage periods.

The use of frozen blood is recommended since freezing and thawing enhance blood cells lysis.

### AMPLIFICATION

Keep all the reagents for the amplification on ice during the execution of the whole procedure. Avoid repeated freeze-thaw cycles of each PCR mix

#### Master Mix preparation (two separated reactions for each sample)

Reagents	Volume/sample
<b>PCR Mix 1</b>	45,75 µl
DNA polymerase (5U/µl), orange cap	0,25 µl

Reagents	Volume/sample
<b>PCR Mix 2</b>	45,75 µl
DNA polymerase (5U/µl), orange cap	0,25 µl

- Prepare each Master Mix for the number of samples extracted + 1 more volume (if “n” ≤ 10) or + 2 more volumes (if “n” > 10).
- Mix gently and dispense **46 µl** of Master Mix in the previously marked 0,2 ml test tubes. Discard the remaining Master Mix.
- Add **4 µl** of extracted DNA to each tube and mix pipetting up and down (**DNA concentration to be used should be 30-40 ng/µl**).

#### Amplification program

- Preheat the thermal cycler lid.
- Put reaction tubes in the thermal cycler and run the following amplification program:

Temperature	Time	Cycles
95°C	2 min.	1
95°C	30 sec	35
60°C	30 sec	
72°C	45 sec	
72°C	5 min	1
10°C	∞	-

#### for C1000 (BioRad) only

Temperature	Time	Cycles
95°C	2 min.	1
95°C	30 sec	38
<b>58°C</b>	30 sec	
72°C	45 sec	
72°C	5 min	1
10°C	∞	-

Keep PCR products at + 2/+8°C if detection is done in the same day, or store them at -25/-15°C.

**Optional:** analyze amplification products by 2% agarose gel containing an intercalating DNA agent.

Mix1	Amplicon Length (bp)
Fattore V R506Q	146
Fattore V H1299R	174
Fattore II G20210A	207
MTHFR C677T	116
MTHFR A1298C	175
CBS 844ins68	242/310
PAI-1	177
GPIIIa	153
ACE I/D	216/494
ATR-1	280
FGB	208
AGT	175
Mix 2	Amplicon Length (bp)
ApoE	278
Fattore XIII	167

## AUTOMATIC DETECTION

Detection can be performed by using automatic instrument (ProfiBlot T48, Dynablot or similar instrument): select the **COAGUL 5** program. In this case **30 µl of mix1 PCR product + 30 µl of mix 2 PCR product + 60 µl DNAT** and 1,5 ml of each reagent should be used for each sample.

## MANUAL DETECTION

- Adjust the water level of the shaking bath to approx. 2/3 of the height of the Plastic Tray.
- Heat the water bath to **45°C** and check the temperature ( $45^{\circ}\text{C} \pm 0,5^{\circ}\text{C}$ ) using a calibrated thermometer.
- Pre-warm the Hybridization and Wash Solution A to **45°C**; all crystals should be completely dissolved.
- Allow DNAT, Conjugate, Wash Solution B and Color developer to reach room temperature.
- Remove one strip for each sample using clean tweezers (touch strips with gloves only). Label strips with a pencil (do not use ballpoint pen, markers, etc).

**Warning:** do not allow the strips to dry during the whole procedure.

## Hybridization (45°C)

- Pipette **40 µl** of **DNAT** into each lane of the Plastic Tray  
**Warning:** do not use DNAT if is not blue.
- For each sample add **20 µl** of **mix 1 amplification product + 20 µl** of **mix 2 amplification product** to the corresponding drop of DNAT. Mix thoroughly with a pipette (solution will remain blue).
- Incubate **5 min** at room temperature.
- Add **1 ml** of **Hybridization** (pre-warmed to 45°C) to each sample and mix gently: blue color disappears.
- Place strips into the respective lanes with marker lines facing up. Submerge completely.
- Incubate for **30 min** at **45°C** in the shaking water bath (about 50 rpm). Keep the water bath closed to avoid temperature variations.

### Wash solution A (45°C)

- After hybridization step remove the tray from the water bath; angle the tray upwards and aspirate the solution from each lane using a pipette or a vacuum aspiration apparatus, without touching the strip surface to avoid its damage.
- Add **1 ml** of **Wash solution A** (pre-warmed at 45°C) to each lane and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution A** (pre-warmed at 45°C).
- Incubate **10 min** at **45°C** in the shaking water bath (about 50 rpm).  
**Warning:** put Wash solution A back in the water bath during incubation
- Remove the solution and add **1 ml** of **Wash solution A** (pre-warmed at 45°C).
- Incubate **10 min** at **45°C** in the shaking water bath (about 50 rpm).
- Remove the solution from each lane.

### Color development (room temperature)

All the next steps have to be performed at room temperature with shaking. Since water takes too much time to reach room temperature, place a tray (for example of polystyrene) on the lid of the water bath to isolate the Plastic Trays from the hot water and use the agitation function of the water bath only.

- Add **1 ml** of **Conjugate**.
- Incubate **20 minutes** while shaking (about 50 rpm) at room temperature.
- Remove the solution, add **1 ml** of **Wash solution B** and rinse briefly.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Wash solution B**.
- Incubate for **5 min** while shaking (about 50 rpm) at room temperature.
- Remove the solution and add **1 ml** of **Color developer**.
- Incubate for **20 min** in the dark, while shaking (about 50 rpm) at room temperature.
- Remove color developer and wash strips with distilled water.
- Let strips dry in the dark on absorbent paper before proceeding with the interpretation of results.

## RESULTS INTERPRETATION

The strip scheme below shows the position of the different probes on the strip.


A line is considered positive when a purple band appears at the end of the detection.

To identify the genotype of a sample, analyze each strip with the help of the transparent interpretation card. In order to correctly align strip and transparent interpretation card, colored marker lines on the top and on the bottom have to be used.

### Warning:

- the color intensity could be different between different lines on the same strip but this doesn't affect the interpretation of the results.
- The detection control line (the upper one) indicates the correct reaction of both conjugate and color developer: it **always** has to stain positive.

## Strip scheme

For each mutation one of the following patterns should be obtained:

- Only wild type line Normal genotype (**homozygous wild type**)
- Wild Type and Mutant lines **Heterozygous** genotype
- Only mutant line **Homozygous mutant** genotype

### Warning:

- For each polymorphism the two probes on the strip correspond respectively to the most (up) and the less (down) common variant;
- The CTR line of CBS is always present, therefore if a positive signal is obtained with the INS line it's not possible to determine whether the sample is heterozygous or homozygous for the CBS insertion.

⇒ - the table below shows the different definitions present in the literature for the GPIIIa target:


Gene	Variants	
ITGB3	GPIIIa var T	GPIIIa var C
	HPA-1 var 1a	HPA-1 var 1b
	PI(A1)	PI(A2)

- The following scheme allows to assign ApoE genotypes:

Nucleotide	Lines on the strip	e3/e3	e3/e4	e2/e3	e2/e4	e2/e2	e4/e4
<b>334T</b>	19) 112T	X	X	X	X	X	
<b>334C</b>	20) 112C		X		X		X
<b>472T</b>	21) 158T			X	X	X	
<b>472C</b>	22) 158C	X	X	X	X		X

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## **TROUBLESHOOTING**

### **FALSE NEGATIVE OR LOW SIGNALS**

- A wrong amount of amplification product has been used. PCR products concentration could be too low due to an inefficient amplification.
- Low signals (except Control Line) and false negative results could be caused by too high temperature during hybridization and/or Wash Solution A steps.

#### **Check carefully water bath temperature**

- Hybridization and Wash Solution A have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.

### **FALSE POSITIVE SIGNALS**

- Non-specific signals could be caused by too low temperature during hybridization and/or Wash Solution A step

#### **Check carefully water bath temperature**

- Hybridization and Wash Solution A have not been well pre-warmed or mixed. Discard remaining reagents and repeat tests with new solutions.
- Water level of the water bath is too low: be sure that level reaches approximately 2/3 of the height of the Plastic Tray
- Check carefully the water bath lid during the Hybridization or Wash Solution A, keep it closed.
- Keep Wash Solution A in the water bath during the incubation steps at **45°C**

### **SINGLE LINE UNHOMOGENEOUS COLORATION**

Shaking rate during detection steps is too low. Check that strips are totally submerged.

## **PERFORMANCE CHARACTERISTICS**

### **Diagnostic Sensitivity**

The diagnostic sensitivity of the kit CVD14 was assessed by analyzing 43 DNA samples positive for the mutations and polymorphisms investigated by the kit.

All the samples were correctly genotyped, resulting in a diagnostic sensitivity of 100%.

### **Analytical Sensitivity**

The analytical sensitivity of the kit CVD 14, expressed as the minimum quantity of the target that can be detected, is equal to 5 ng/μl DNA.

### **Diagnostic Specificity**

The diagnostic specificity of the kit CVD 14 was determined by analyzing 43 samples. All samples were correctly amplified and genotyped. In agreement with these data, the diagnostic specificity is 100%.

### **Reproducibility**

#### **INTRA-ASSAY**

The intra-assay reproducibility was assessed by analyzing 1 sample in 28 replicates in the same run and the same samples in 25 replicates in a second run. Two different operators tested the replicates using two batch lots.

#### **INTER-ASSAY**

The inter-assay reproducibility was tested using at least 3 samples in 2 replicates. Two different operators tested several replicates in different runs and days, using same lot of reagents. The strips development was made in automatic and manual.

Tests were performed as described in the instruction for use and the reproducibility was 100%.

## WARNING

As every system based on amplification and detection of nucleic acids, it is possible to have unexpected results in the presence of unknown genetic mutation in the primer and/or probe specific target regions.